

5L Medical Oxygen Concentrator

D02-5E



These models can produce average 93% purity oxygen to meet 5liter flow requirements of hospital, clinic, home use, nursing centre. It can stably work in 24 hours a day, always maintain oxygen purity at average 93% at 5L/min flow. And in future, it only needs few maintenance of the main machine among its whole lifetime.

5L Medical Oxygen Concentrator

D02-5E



Technical data	
Oxygen flow	1-5 LPM
Output pressure	50-70 kPa
Power supply	AC 230V, 50Hz or 115V, 60Hz
Input power	300VA (AVER)
Oxygen concentration	1-5 LPM 93%±3%
Operation Environment	Temperature 5° C ~40° C
	Humidity≤80%
	Atmospheric pressure:86kPa~106 kPa
OCSI indicator	>82% green; <82% red
N.W.	14kg
Noise	42dB(A)
Size	38cm×26cm×58 cm(L×W×H)
Operation regular	Time cycle/pressure rotation
Transportation/storage condition	Temperature-20° C ~55° C
	Humidity≤90%
	Atmospheric pressure:86kPa~106 kPa
Equipment class&Type	Class II B



Email: info@dynmed.cn

Floor 2 of No .2 Workshop, No. 66, Shunxing Road, Lujia Town,
Kunshan, Suzhou City, Jiangsu Province, 215331, China